

Pro-MicroScan

Microscope Digital Camera

BASIC USER MANUAL
Version 3.0

USB 2.0 CAMERA



The lightning flash with an arrowhead symbol within an equilateral triangle is intended to alert the user of the presence of uninsulated "dangerous voltage" within the product's enclosure that may be of sufficient magnitude to constitute a risk of electric shock to persons.



The exclamation point within an equilateral triangle is intended to alert the user of the presence of important operating and maintenance (servicing) instructions in the literature accompanying the appliance.

WARNING!

TO AVOID THE RISK OF FIRE OR ELECTRICAL SHOCK, NEVER EXPOSE THIS PRODUCT TO WATER OR OPERATE IN A HIGH HUMIDITY ENVIRONMENT.

Keep camera away from the following:

- High temperatures and excessive humidity
- Direct sunlight or other heat sources
- Dust
- Extreme vibration

If the camera lens becomes dirty, first try blowing dust particles from it. Next use a lint free cloth or proper lens cleaning cloth to wipe the lens gently. Only use approved lens cleaning fluids. Never touch the lens with your fingers! Do not allow the camera to come into contact with water or any other liquid. If water or a foreign substance should get in the camera, unplug the camera immediately and do not attempt to plug in until after a long period of time to allow for drying. This does not necessarily mean it will work again! Do not open the camera case or attempt your own repairs. Internal components may create risk of electrical shock or fire. There are no user serviceable components inside. Clean the outside of the camera by wiping with a clean, dry cloth. Never use harsh or abrasive cleaners or organic solvents on the camera or any of its parts. If you notice smoke or a strange smell coming from the camera, immediately unplug the camera and consult an authorized dealer.

Trademark and Copyright information

© COPYRIGHT 2008, OPLENIC

All rights are reserved. No part of this document may be photocopied, reproduced, or translated to another language without the prior written consent of OPLENIC. The information contained in this document is subject to change without notice. OPLENIC shall not be liable for errors contained herein or for consequential damages in connection with the furnishing, performance or use of this material.

Microsoft[®] and Windows[®] are U.S. registered trademarks of Microsoft Corporation.

Welcome

Thank you for buying OPLENIC ***Pro-MicroScan*** high resolution USB2.0 camera. We suggest you spend a few minutes carefully reading this user manual before installing and using your camera.

Please keep this user manual in a safe place for future reference.

Warnings and Precautions

Precautions when using your *Pro-MicroScan*



Keep camera away from the following:

- High temperatures and excessive humidity
- Direct sunlight or other heat sources
- Dust
- Extreme vibration

If the imaging sensor lens becomes dirty, first try blowing dust particles from it. Next use a lint free cloth or proper lens cleaning cloth to wipe the lens gently. Only use approved lens cleaning fluids. **Never touch the lens with your fingers!**

Do not allow the camera to come into contact with water or any other liquid. If water or a foreign substance should get in the camera, unplug the camera immediately and do not attempt to plug in until after a long period of time to allow for drying. This does not necessarily mean it will work again! Do not open the camera case or attempt your own repairs. Internal components may create risk of electrical shock or fire. There are no user serviceable components inside.

Clean the outside of the camera by wiping with a clean, dry cloth. Never use harsh or abrasive cleaners or organic solvents on the camera or any of its parts.

If you notice smoke or a strange smell coming from the camera, immediately unplug the camera and consult an authorized dealer.

System Requirements

PC running either Windows[®] XP, Vista, 2008

- ◆ DirectX 9.0 or later installed
- ◆ 1.8 GHz minimum processor
- ◆ USB 2.0 Port
- ◆ CD-ROM drive
- ◆ 512 MB RAM(Recommend 1GB or more)
- ◆ 200 MB Free HDD

Performance is dependent on the PC's graphics card.

USB 2.0 is only supported on Microsoft Windows[®] XP, 2003, Vista and 2008. Do not try to use any other Microsoft[®] operating system. will NOT work on a Mac Computer.


Windows XP Professional must have:

Service Pack 2 (or later) installed in order for the camera to function correctly. To check if you have USB2.0 installed on your computer goes to:

Start > Settings > Control Panel > System... or 'right-click' on "My Computer" and select "properties" from the context menu...


Click on the "Hardware" tab and then onto the button "Device Manager";

Within Universal Serial Bus controllers (USB controllers) you will see a list of USB devices installed on your machine. In this example the red arrow is indicating the USB2.0 device. USB2.0 controllers may be labeled as USB2.0 controllers, however they may also be labeled as "Enhanced" or "Standard Enhanced" controllers. The so-called "Standard" or "Standard Universal" controllers are for the slower USB1.x system.



Package Contents

If you find any contents missing or damaged, contact the place of purchase



- 1 PCS EYEPiece CAMERA
- 1 PCS DUST COVER
- 1 PCS MINI USB2.0 CABLE
- 1 PCS SOFTWARE CD (OPTION)

*Software can be download from website DEM page for free.


Installing the Software

(If you do not have the software CD, skip this step.)

1. Insert CD into CD-ROM drive
2. Install image capture software (After installation the image software icon will created on the desktop)



To Connect Camera

Connect the USB cable into a free USB slot on the computer (Requires **USB 2.0** interface). 

“New Hardware Found” displays at the bottom right of your PC screen. The camera is recognized.


NOTE:



1. Plug the USB connector to board USB2.0 port directly on the back of the PC is recommend,
2. Do NOT use anything to extend the USB cable.

Using the camera

There are 2 ways to get video from the camera.

A. Using Microsoft Windows Viewer

1. Open “My Computer”, you will find the “USB Video Device” with  icon.



2. Double Click the “USB Video Device” icon to start video.

3. Capture and delete...

B Using video image software

1. Start the MiniSee

2. Start video


MiniSee--Click on the camera name below “Video” button.

3. Connect the camera to the microscope and turn on the light.


(Please make sure your light system has the option of a BLUE LIGHT FILTER. For the CCD to view properly, you would need light color temperatures around 5600K-6000K.)

4. Setting the Video source property

5. Capture the image Press  to capture image



6. Capture the Video



Trouble Shooting

WARNING:

Please do not attempt to open the camera case. There are no user serviceable parts inside.

USB Device Conflicts

When there are a number of USB devices installed at the same time, unexpected errors can arise. The cause can be insufficient available USB bandwidth, or possibly too much combined current being drawn by USB devices. To test for these conditions you should uninstall and unplug as many other USB devices as possible. This should help to rule out the above, and eliminate possible conflicts with other cameras. If you are using an extended USB cable, please make sure the cable's bandwidth is good enough to support High-Speed USB 2.0 standard. Once the camera is not working well, please remove it.

Image Preview/Capture Problems

If you have installed camera and you get a black screen, Check the setting of camera, make sure preview is select on

Frame rate slow & color noise on the screen

Video Mode: use smaller size video mode will get faster frame rate video

Integration time / luminous intensity setting: longer Integration time and brighter illumination will make the frame rate faster.

If you have any other problems, please visit our website to get help information

<http://www.oplenic.com/support.asp>

<http://www.oplenic.com/faq.asp>

Technical Data Sheet

Pro-MicroScan : SkopKam digital camera FOR microscope

Model	DEM35	DEM130	DEM200	DEM300	DEM500
Part Number	#5820	#5821	#5822	#5823	#5825
Sensitivity (v/lux-sec @550nm)	2.0	1.0	2.0	1.0	1.0
Max Resolution (Hardware)	640X480 0.3M pixels	1280X1024 1.3M pixels	1600X1200 2M pixels	2048X1536 3M pixels	2592 X 1944 5M pixels
Full Res. Frame Rate (Depend on the PC)	30frames/sec@ 640X480	15frames/sec@ 1.3M pixels	15frames/sec@ 1600X1200	3frames/sec@ 2048X1536	3frames/sec@ 2592 X 1944
Dynamic Range	60dB	68dB	71dB	75dB	75dB
Wave Length(nm)	400-650nm (IR-Cut added)				
Exposure	ERS (Electronic Rolling Snap)				
White Balance	Automatic / Advanced Manual				
Auto-exposure	Automatic / Advanced Manual				
Image Output	USB2.0, 480Mb/s				
Power Supply	USB2.0 DC 5V / 250mA				
Programmable Control	Image size, exposure time, brightness, gain, Gamma Correction, color,, contrast, Saturation,				
Working Temperature	-30°C ~ 70°C				
USB Cable	Mini USB cable, length 1.5m				
Accessories Included	2 adaptors (Φ23.2 to Φ30, Φ23.2 to Φ30.50, one each) available for option at extra cost				
Output data	RGB24, VFW,Direct Show and Twain interface				
Software	Image software: MiniSee				
Operating System	Windows 2000 / Me / XP(SP2/SP3) / 2003 / Vista / 2008 32Bit system				

Warranty

OPLENIC Three Year Limited Warranty

OPLENIC warrants its products to be free from the date of purchase this product for three year. OPLENIC will repair or replace such product or part thereof which, upon inspection by OPLENIC, is found to be defective in materials or workmanship. As a condition to the obligation of OPLENIC to repair or replace such product, the product must be returned to OPLENIC together with proof-of-purchase satisfactory to OPLENIC.



All returns must be accompanied by a written statement setting forth the name, address and daytime telephone number of the owner, together with a brief description of any claimed defects. Parts or product for which replacement is made shall become the property of OPLENIC.

The customer shall be responsible for all costs of transportation and insurance, both to and from the factory of OPLENIC, and shall be required to prepay such costs.

OPLENIC shall use reasonable effort to repair or replace any accessory covered by this warranty within thirty days of receipt. In the event repair or replacement shall require more than thirty days, OPLENIC shall notify the customer accordingly. OPLENIC reserves the right to replace any product which has been discontinued from its product line with a new product of comparable value and function.

License Agreement

The software and accompanying documentation are protected by the International copyright laws. Any use of this software in violation of copyright laws or the terms of this agreement will be prosecuted to the best of our ability.

OPLENIC authorizes you to make archival copies of this software for the sole purpose of back-up and protecting your investment from loss. Under no circumstances may you copy this software or documentation for the purposes of distribution to others. You may not distribute, modify, adapt, translate, reverse engineer, recompile, disassemble or otherwise attempt to discover the source code of the Software. Any information supplied by OPLENIC or obtained by you, as permitted hereunder, may only be used by and may not be disclosed to any third party or used to create any software which is substantially similar to the Software.